

Indication

The Disposable Powered Endoscopic Linear Stapler & Reloads are provided for multiple open or minimally invasive surgical procedures for anastomotic creation and closure of stump or incision.

Product Parameters

1、Device Body Parameters


Model	Length of arm(mm)	Articulation Angle(°)	Model	Length of arm(mm)	Articulation Angle(°)
CNPBS45	190±5	45°±10°	CNPPS45	190±5	45°±10°
CNPBS60	190±5	45°±10°	CNPPS60	190±5	45°±10°
CNPBM45	250±5	45°±10°	CNPPM45	250±5	45°±10°
CNPBM60	250±5	45°±10°	CNPPM60	250±5	45°±10°
CNPBL45	350±5	45°±10°	CNPPL45	350±5	45°±10°
CNPBL60	350±5	45°±10°	CNPPL60	350±5	45°±10°

2、Reloads Parameters

Model	Staple Height	Number of Staples	Suture Length (mm)	Reload Color	Model	Staple Height	Number of Staples	Suture Length (mm)	Reload Color
CNERSA45	2.0±0.2	70	49.3±2	Gray	CNERPA45	2.0±0.2	70	49.3±2	Gray
CNERSW45	2.5±0.2	70	49.3±2	White	CNERPW45	2.5±0.2	70	49.3±2	White
CNERSB45	3.5±0.2	70	49.3±2	Blue	CNERPB45	3.5±0.2	70	49.3±2	Blue
CNERSY45	3.8±0.1	70	49.3±2	Golden	CNERPY45	3.8±0.1	70	49.3±2	Golden
CNERSG45	4.1±0.1	70	49.3±2	Green	CNERPG45	4.1±0.1	70	49.3±2	Green
CNERSH45	4.3±0.1	70	49.3±2	Black	CNERPH45	4.3±0.1	70	49.3±2	Black
CNERSA60	2.0±0.2	88	61.3±2	Gray	CNERPA60	2.0±0.2	88	61.3±2	Gray
CNERSW60	2.5±0.2	88	61.3±2	White	CNERPW60	2.5±0.2	88	61.3±2	White
CNERSB60	3.5±0.2	88	61.3±2	Blue	CNERPB60	3.5±0.2	88	61.3±2	Blue
CNERSY60	3.8±0.1	88	61.3±2	Golden	CNERPY60	3.8±0.1	88	61.3±2	Golden
CNERSG60	4.1±0.1	88	61.3±2	Green	CNERPG60	4.1±0.1	88	61.3±2	Green
CNERSH60	4.3±0.1	88	61.3±2	Black	CNERPH60	4.3±0.1	88	61.3±2	Black
CNERPM45	3.8,4.1,4.3	70	49.3±2	Purple	CNERPL60	2.5,3.0,3.5	88	61.3±2	Brown
CNERPL45	2.5,3.0,3.5	70	49.3±2	Brown	CNERPM60	3.8,4.1,4.3	88	61.3±2	Purple



Disposable Powered Articulating Endoscopic Linear Cutter Stapler



Jiangsu CHANNEL Medical Device Co., Ltd.

Tel: +86-512-36857161
Web: www.channel-mdi.com
Email: marketing@channel-mdi.com
Add: Fourth Floor, G3, 999#, Qujia Road, Qiandeng Town, Kunshan, Jiangsu, PEOPLE’S REPUBLIC OF CHINA

